Technical Data Sheet



THE SOFT MATERIAL FOR GINGIVAL MASKS **FEATURES**

GL 4000 is a gingiva-like material designed for the production of realistic models of soft tissues and gingival masks. Its pink color and its consistence remind the natural gingiva ones.

This material is developed and produced in-house by DWS.

- Smooth Surface
- High Resolution
- High Accuracy
- Soft Material
- Pink Colour

UV post treatment after the print and the cleaning is not required.

TECHNICAL CHARACTERISTICS OF THE LIQUID MATERIAL

Environmental Values for Use	22°C - 27°C - max, RH 40% - 60%
Appearance / Colour	Liquid / Pink
Density	1,00 g/cm ³
Viscosity	1400 ~ 2100 mPa•s at 25°C

TECHNICAL CHARACTERISTICS OF THE RESIN AFTER UV CURING

Elongation at Break (%)	25 ~ 65
Surface Hardness (Shore D)	20 ~ 24
Application / Use	Soft Tissue Models, Gingival Masks

Technical specification subject to change without notice.



